

4dentPRO
white label solutions

About us 1-2

Prevention 3-4

Restoratives 5-7

Impression Material 9-13

Bite Registration 15-16

Soft Relining 17-18





Perfect fit. Innovative. Versatile.

- ✓ In-house research and product development
- ✓ Quality testing and product release according to the latest legal requirements
- ✓ Expert support in building your brand
- ✓ Guaranteed long-term product availability
- ✓ Organization and management of your product packages
- ✓ Handling of worldwide transportation
- ✓ Support in order planning to adhere with delivery deadlines

Just three steps to your finished, desired product!

1.

In the beginning there is your idea!

It is our reference and basis for in-depth consulting & conception, resulting in a detailed basis for planning your desired end product. A competent team with comprehensive technical expertise and many years of expert experience is at your disposal to realize your special requirements.

We find the ideal product solution

2.

Individual approval sample for testing? No problem!

On request, you will receive a sample of your product selection in advance with the agreed specifications for internal testing and evaluation, so we really can be sure to meet your requirements. For 4dentPro, sampling is an essential precursor for product finalization to meet customer needs.

You have the chance to check an individual approval sample in advance

3.

Lean back, we'll take care of the rest!

Following extensive evaluation of our samples, you will receive the final product specifications along with all contract documents. Besides fast and uncomplicated product finalization, our team of experts handles all formalities for marketing authorization in close coordination with you and supports you competently in handling your delivery and organizing transportation.

Product finalization with optional add-on services.



Prevention

FLUORIDE VARNISH – THE NEW DIMENSION IN PREVENTION

- ✓ Use as part of professional tooth cleaning
- ✓ Desensitization
- ✓ Easy handling
- ✓ Various flavors available

Prevention

Fluoride Varnish



Overview

Fluoride Varnish adheres immediately and securely to the tooth surface for several hours → thanks to the unique dimethicone gel base.



User benefits

- No prior drying of the teeth necessary
- Can be applied in any direction and reaches all areas
- A single dose is sufficient for treating all teeth in both dental arches



Patient benefits

- Pain relief after treating hypersensitivities
- Protective layer for the tooth surface
- No burning and no bitter aftertaste >> without ethyl alcohol
- Non-allergenic formula
- Different flavors available



Minimum order quantity:

Blister	0,5ml	93.000 Blister
---------	-------	----------------



Restoratives

SOLUTIONS FOR CORE BUILD-UP, ROOT POST FIXATION, TEMPORARY CROWNS AND BRIDGES, LUTING COMPOSITE

- ✓ Free of bisphenol A
- ✓ Easy handling
- ✓ Stable and fracture resistant
- ✓ Different container types and accessories



Restoratives



Temporary Crown and Bridge Material



Overview

Temporary crown and bridge material based on a multifunctional acrylic composite, suitable for the fabrication of temporary crowns, partial crowns, bridges, inlays, onlays and veneers including long-term provisional restorations.

different cartridge system available:

10:1 cartridge system shades: A1, A2, A3, A3.5, B1, BL

1:1 cartridge system shades: A1, A2, A3, BL



Benefits

- Very stable and fracture-resistant
- Only minimal dust created during milling
- Minimal inhibition layer after curing
- High fluorescence and translucency



Minimum order quantity:

cartridge	1 x 50 ml	300 units
syringe	1 x 10 ml	1,500 units

Core Build-up Material



Overview

Core build-ups are compression-resistant and form a reliable monoblock with the root post. Indicated for building up cores and cementing root posts. Available in white, dentin and blue.



Benefits

- Can be applied directly into the cavity
- Easy and safe handling
- Flows into the root canal well
- Dual curing ensuring strong and reliable core build-ups even in areas inaccessible to light



Minimum order quantity:

syringe	1 x 2.5 ml	1,200 units
syringe	1 x 5 ml	500 units
cartridge	1 x 25 ml	100 units

Restoratives

Dual-curing luting composite



Overview

Dual-curing, adhesive luting composite suitable for cementing all indirect restorations and compatible with all dental materials (metal, oxide ceramic, etc.). Luting composite can be used to cement conventional restorations, such as crowns and bridges, as well as challenging and highly esthetic restorations such as adhesive bridges and veneers. Luting composite can also be used to create core build-ups.



Benefits

- Two-in-one: luting and core build-up
- Precise and reliable results
- High color stability
- Easy and safe handling



Minimum order quantity:

syringe	1 x 2.5 ml	1,200 units
syringe	1 x 5 ml	500 units
Tooth Primer	1 x 4 ml	
Restorative Primer	1 x 4 ml	



Impression material

FOR EVERY IMPRESSION TECHNIQUE THE PERFECT PRODUCT

- ✓ Different setting times
- ✓ Different viscosities
- ✓ Simple handling
- ✓ Wide variety of colors and flavors
- ✓ Different container types and accessories



A-Silicone Putty (regular/fast/soft)



Overview

A kneadable and addition-curing preliminary impression material based on vinyl polysiloxane for taking preliminary impressions with different degrees of final hardness and setting characteristics.



Benefits

- Various setting times and final hardnesses
- Easy handling
- Top elastomeric properties to ensure easy and reliable removal from the mouth
- Neutral odor and taste



Minimum order quantity:

Jars	2 x 450 ml	600 units
------	------------	-----------

A-Silicone Jumbo Cartridge Putty (fast/soft)



Overview

A kneadable and addition-curing precision impression material based on vinyl polysiloxane for use when fabricating inlays, onlays and precision-fit crowns and bridges.



Benefits

- Genuine putty in a 5:1 cartridge
- Easy handling and hygienic mixing
- Reduced material wastage
- Neutral odor and taste



Minimum order quantity:

Jumbo cartridge	1 x 380 ml	1,000 units
-----------------	------------	-------------

A-Silicone Heavy Body (fast/soft)



Overview

A high viscosity and addition-curing precision impression material based on vinyl polysiloxane with different degrees of final hardness and setting characteristics.



Benefits

- User- and patient-friendly setting times
- High tear-resistance
- Top elastomeric properties
- Neutral odor and taste



Minimum order quantity:

Jumbo cartridge	1 x 380 ml	1,000 units
Cartridge	1 x 50 ml	5,500 units

A-Silicone Light/X-Light Body (regular/fast)



Overview

A low viscosity, addition-curing and elastomeric precision impression material based on vinyl polysiloxane. Ideal in combination with a high viscosity material for the one-step and two-step putty-wash techniques.



Benefits

- User- and patient-friendly setting times
- Good thixotropic properties
- Flows precisely and accurately even under extreme conditions
- Neutral odor and taste



Minimum order quantity:

Cartridge	1 x 50 ml	5,500 units
-----------	-----------	-------------

A-Silicone Medium Body (regular/fast)



Overview

A medium viscosity and elastomeric monophasic impression material based on A-silicone for taking precision impressions and impressions for fixed or removable prosthetic restorations and in implantology.



Benefits

- Ideal intraoral setting time in addition to user-friendly total working times
- High thixotropy
- Optimum flow properties
- Neutral odor and taste



Minimum order quantity:

Jumbo cartridge	1 x 380 ml	1,000 units
Cartridge	1 x 50 ml	5,500 units

A-Silicone Lab Putty



Overview

A kneadable and addition-curing overcast material with high final hardness, used as overcast and bite index material and also ideal for creating dental-technical constructions.



Benefits

- High final hardness: precise fixation on the overcast
- Dimensionally accurate elastic recovery
- Cleanly and easily dispensed
- Soft kneading without sticking



Minimum order quantity:

Bucket	2 x 1 kg	500 units
--------	----------	-----------

Impression material

Alginate Alternative



Overview

A medium viscosity A-Silicone with alginate-like consistency for taking impressions of opposing jaws to fabricate temporary crowns and bridges as well as splints.



Benefits

- Optimum elastic recovery
- High dimensional stability
- High tear-resistance
- Neutral odor and taste



Minimum order quantity:

Jumbo cartridge	1 x 380 ml	1,000 units
Cartridge	1 x 50 ml	6,500 units

Polyether-Alternative (medium/heavy/light)



Overview

Elastomeric impression material based on a combination of vinyl polysiloxane and polyether available in different viscosities and setting characteristics for taking precision impressions.



Benefits

- User- and patient-friendly setting times
- Available in high, medium and low viscosity
- Easy handling
- Neutral odor and taste



Minimum order quantity:

Jumbo cartridge	1 x 380 ml	1,000 units
Cartridge	1 x 50 ml	5,500 units



Bite registration

RIGHT PRODUCT FOR ALL REQUIREMENTS

- ✓ Bite registration for precise capture of the occlusion conditions
- ✓ Precise results with easy handling
- ✓ Different setting times and final hardness
- ✓ Wide variety of colors and flavors
- ✓ Different container types and accessories



Bite Registration



Bite Registration/Extra Hard (regular/fast)



Overview

An addition-curing bite registration material based on vinyl polysiloxane and with high final hardness.



Benefits

- Ideal elasticity
- User- and patient-friendly setting times
- High thixotropy
- Precise results



Minimum order quantity:

Cartridge	1 x 50 ml	5,500 units
-----------	-----------	-------------



Bite Registration Clear (regular/fast)



Overview

An addition-curing, transparent bite registration material based on vinyl polysiloxane.



Benefits

- Guaranteed transparency for optimal aesthetics
- Precise results
- Highly convenient
- Neutral odor and taste



Minimum order quantity:

Cartridge	1 x 50 ml	5,500 units
-----------	-----------	-------------



Soft Relining



Soft Relining

FOR DIRECT AND INDIRECT USE

- ✓ A-silicone-based Soft Relining material
- ✓ Permanently soft over long periods of wearing
- ✓ Easy handling
- ✓ Neutral odor and taste
- ✓ Different container types and accessories



Permanently Soft Relining Material



Overview

Permanently soft relining material based on vinyl polysiloxane for direct and indirect application.



Benefits

- Excellent bond-strength on the denture
- Very smooth and hydrophobic silicone surface
- Permanently elastic
- High tear-resistance



Minimum order quantity:

Soft Relining	1 x 50 ml	500 units
Silicone Sealant	1 x 50 ml	500 units
Starter Set, finished sales unit	<ul style="list-style-type: none"> • soft relining material, • silicone sealant, • adhesive, • mixing tips 6.5 mm, • mixing tips 3.2 mm, • brush holder, • single use brushes, • accessories 	500 sets

